



cran en cuivre plaqué tain 6 curs isol en PFA, cble / fil de capteur PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

cran en cuivre plaqu tain 6 curs isol en PFA, cble / fil de capteur PRT

Temperature leve jusqu 260C, construction en tflon PFA extrud. Couramment utilis pour les thermomtres Pt100, galement adapt dautres types tels que Pt130, Pt500 et Pt1000

Caractristiques techniques**Specifications**

General Description	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores with nylon cord, tin plated copper braided screen with PFA insulated outer jacket
Cable Length	10m, 25m, 50m, 100m
Cable Type	PFA Teflon® construction
Core / Strands	7/0.2mm
Insulation Material	PFA Insulated
Max. Temperature	+260°C
Min. Temperature	-75°C
Number of Wires	6
Standards Met	Cores colour coded Red & White in accordance with IEC-751, black outer jacket