

## Thermopockets for $\nabla f \hat{\text{A}}^{\text{C}}$ s solides



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Les thermopockets en acier inoxydable de différentes tailles sont conçus pour une utilisation avec des sondes de température de style industriel afin de protéger la gaine des supports corrosifs et de faciliter le remplacement des sondes sans perturber le processus.

Du stock, les connexions de processus sont BSPP de 1/2 pouce en standard, un raccord de compression BSPP de 1/8 pouce pour convenir au diamètre de la sonde est nécessaire pour connecter le capteur au thermopocket.

## Specifications

<b>General Description</b>	From stock, process connections are 1/2 inch BSPP as standard, a 1/8 inch BSPP compression fitting to suit diameter of probe is required to connect sensor to thermopocket.
<b>Material</b>	Stainless steel
<b>Dimensions</b>	303SS, 12mm Stem
<b>Type &amp; Stem Diameter</b>	303SS, 12mm Stem
<b>Stem Length</b>	50mm, 100mm, 150mm
<b>Probe Diameter</b>	50mm, 100mm, 150mm stem length