



Thermopockets fabriquÃ©s



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Fabricated Thermopockets in varying sizes for use with Industrial Style Temperature Probes.
Material	Stainless steel
Length	316SS with 7/16" or 8mm stem
Type & Stem Diameter	316SS with 7/16" or 8mm stem
Stem Length	50mm, 100mm, 150mm, 350mm, 550mm
Probe Diameter	100mm, 150mm, 350mm, 550mm