




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple jonction expose IEC double torsade isol en PTFE avec prise miniature - Type K - Longueur 1M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple de type K, J ou T. Rponse rapide, jonction expose soude disponible dans la spcification de cble de classe 1 conformment la norme BS EN 60584-3:2008 avec une terminaison de fiche miniature.

Idal pour les applications de test et de developpement.

Spfications

Specifications

Product Code	XE-3505-001
General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	K IEC
Cable Length	1M
Cable Type	Teflon® PTFE insulated twin-twisted lead
Core / Strands	1/0.2mm
Termination Type	Miniature flat pin plug colour coded in accordance IEC 584
Max. Temperature	250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3
Connector Type	Plug
Connector Size	Miniature
Series	Exposed Junctions