

**Thermocouple jonction expose IEC paires plates  
isoles en fibre de verre - Type T avec terminaison  
queues nues - 1/0,315 mm x 2 M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Caractéristiques techniques

**Specifications**

<b>Product Code</b>	XE-3561-001
<b>General Description</b>	Good temperature resistance (not suitable in fluids)
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	T
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre
<b>Core / Strands</b>	1/0.315mm
<b>Termination Type</b>	Bare Tails
<b>Max. Temperature</b>	+350°C (short periods up to 400°C)
<b>Min. Temperature</b>	-60°C
<b>Tolerance</b>	Class 2
<b>Standards Met</b>	IEC-584-3