



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Thermocouple jonction expose IEC isol PFA  
(conducteurs de 0,076 mm) - Type K - Longueur  
de 0,5 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple de type K. Rponse rapide, jonction expose soude disponible avec queues nues. Idal pour les applications de test et de dveloppement.

### Spifications

## Specifications

<b>Product Code</b>	XE-3509-001
<b>General Description</b>	Fast response, welded tip exposed junction
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	K IEC
<b>Cable Length</b>	0.5M
<b>Cable Type</b>	Teflon® PFA insulated wires
<b>Core / Strands</b>	0.003" (0.076mm) conductors
<b>Termination Type</b>	Bare Tails
<b>Max. Temperature</b>	260°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	IEC-584-3
<b>Series</b>	Exposed Junctions