



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple junction expose IEC isol PFA (conducteurs de 0,076 mm) - Type K - Longueur 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple de type K. Rponse rapide, jonction expose soude disponible avec queues nues. Idal pour les applications de test et de dveloppement.

Spfications**Specifications**

Product Code	XE-3511-001
General Description	Fast response, welded tip exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	K IEC
Cable Length	2M
Cable Type	Teflon® PFA insulated wires
Core / Strands	0.003" (0.076mm) conductors
Termination Type	Bare Tails
Max. Temperature	260°C
Min. Temperature	-75°C
Standards Met	IEC-584-3
Series	Exposed Junctions