




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Coussinets de fixation thermique en tissu de verre PTFE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Tampons auto-adhésifs en tissu de verre PTFE pour la localisation des thermocouples pointes soudées

Class -73C 260C

Taille du tampon 13x18mm

La feuille contient 20 tampons

Specifications

Product Code	XE-3599-001
Series	ATTACHMENT PADS