



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Thermocouple isol minral IEC avec bouchon standard - Type K - 1,5 mm x 500 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Vous avez besoin d'atlonner votre capteur de temperature ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility propose le choix entre un talonnage traable en 3 ou 5 points via notre installation dtalonnage interne.

**Veillez noter que si vous ajoutez une option de calibration votre commande, vous serez contact pour les informations de livraison.**

**Caractristiques techniques****Specifications**

<b>Product Code</b>	XE-3702-001
<b>General Description</b>	Highly flexible, sheath can be bent/formed to suit many applications and processes
<b>Sensor Type</b>	Mineral Insulated
<b>Thermocouple Type</b>	K
<b>Element/hot junction</b>	Single element, junction insulated from sheath (offers protection against spurious electrical signals)
<b>Probe Diameter</b>	1.5mm
<b>Probe Length</b>	500mm
<b>Sheath Material</b>	310 Stainless Steel
<b>Termination Type</b>	Standard round pin plug, colour coded in accordance IEC 584-3
<b>Max. Temperature</b>	1100°C
<b>Min. Temperature</b>	-40°C
<b>Standards Met</b>	In accordance with IEC 584-1&2
<b>Connector Type</b>	Plug



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

<b>Connector Size</b>	Standard
<b>Colour</b>	Green
<b>Series</b>	Mineral Insulated