

Thermocouple minral isol de style industriel avec tte terminale KNE - Type K - 6,0 mm de diamtre - 400 mm de longueur



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Vous avez besoin d'étalonner votre capteur de température ?

Traçable selon les normes d'étalonnage ISO17025 (UKAS), Labfacility propose le choix d'un étalonnage traçable en 3 ou 5 points grâce à notre installation d'étalonnage interne ; vous pouvez sélectionner n'importe quels points de température entre -10°C et +200°C.

Veuillez noter que si vous ajoutez une option de calibration à votre commande, vous serez contacté pour les informations de livraison.

Caractéristiques techniques**Specifications**

Product Code	XE-3254-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral insulated Thermocouple fitted with KNE terminal head
Thermocouple Type	K
Thermocouple Junction	Insulated hot junction
Probe Diameter	6.0mm
Probe Length	400mm
Sheath Material	316 stainless steel sheath
Head Type	KNE, Cast aluminium finish, screw type lid with chain
Gland	M16 x 1.5mm Plated brass gland fitted
Max. Temperature	1100°C
Min. Temperature	-40°C
IP Rating	IP67