


**Thermocouple jonction expose isol ANSI PFA  
avec prise miniature - Type K - Taille du brin 1/0,3  
mm - Longueur 2 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

|                            |                                         |
|----------------------------|-----------------------------------------|
| <b>Product Code</b>        | XE-3582-001                             |
| <b>General Description</b> | Fast response, welded exposed junction  |
| <b>Sensor Type</b>         | Exposed Junction                        |
| <b>Thermocouple Type</b>   | K                                       |
| <b>Cable Length</b>        | 2M                                      |
| <b>Cable Type</b>          | Teflon® PFA insulated twin-twisted lead |
| <b>Core / Strands</b>      | 1/0.3mm                                 |
| <b>Termination Type</b>    | Miniature Plug                          |
| <b>Max. Temperature</b>    | 250°C                                   |
| <b>Min. Temperature</b>    | -75°C                                   |
| <b>Standards Met</b>       | ANSI MC96.1                             |
| <b>Series</b>              | Exposed Junctions                       |