



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Thermocouple isolation minrale ANSI - Type K avec fiche miniature - 1,5 mm x 250 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Vous avez besoin de calibrer votre capteur de température ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne.

**Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.**

**Spcifications****Specifications**

<b>Product Code</b>	XE-3001-001
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Hot junction insulated from sheath
<b>Probe Diameter</b>	1.5mm
<b>Probe Length</b>	250mm
<b>Sheath Material</b>	Stainless steel 310 sheath
<b>Termination Type</b>	Miniature Plug
<b>Max. Temperature</b>	1100°C
<b>Standards Met</b>	ANSI