



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple isolation minrale ANSI - Type K avec prise miniature - 3,0 mm x 500 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Vous avez besoin de calibrer votre capteur de température ?

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne.

Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.

Spifications**Specifications**

Product Code	XE-3017-001
Thermocouple Type	K
Thermocouple Junction	Hot junction insulated from sheath
Probe Diameter	3.0mm
Probe Length	500mm
Sheath Material	Stainless steel 310 sheath
Termination Type	Miniature Plug
Max. Temperature	1100°C
Standards Met	ANSI