


Capteur de temprature de lair / Capteur PT100 intrieur



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-3670-001
General Description	PT100 wall mounting indoor temperature sensor for environmental temperature measurement
Sensor Type	For indoor use only
Housing	Moulded case L85 x W85 x H30mm
Max. Temperature	+70°C
Min. Temperature	-50°C
Accuracy	Class B to IEC 751
Number of Wires	Can be connected to any 2, 3 or 4 wire instrument
Features #1	Vented construction
Features #2	Surface or standard electrical conduit mounting
Features #3	Cable entry slots underside of housing
Series	Air Indoor Sensor