



Capteurs de temprature extrieurs / entrepts frigorifiques Pt100 - Type PRT - Dual Element



PT100 OUTDOOR/INDOOR TYPE
TEMPERATURE SENSORS



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Duplex 2 x 4 fils avec deux éléments de détection indépendants, permettant la connexion de deux appareils de mesure ou de contrôle ; alternativement, un capteur peut être utilisé comme capteur primaire et l'autre comme capteur de secours, éliminant ainsi le besoin de retirer et de remplacer la sonde si l'un des capteurs venait tomber en panne - peut être connecté à n'importe quel instrument Pt100 2, 3 ou 4 fils.

Pt100 (100 ohms @ 0°C) selon IEC 751, classe B.

Specifications

Specifications

Product Code	XE-3674-001
General Description	Weatherproof IP67 rated wall mounting Pt100 for accurate outdoor air temperature sensing
Sensor Type	Dual Element: Pt100 (100 Ohms @ 0°C) to IEC 751, Class B
Housing	Painted Die-Cast Aluminium Housing L80 x W74 x H54mm
Element	Pt100, class B to IEC 751
Sheath Material	Perforated 316 stainless steel probe, 12.7mm x 75mm
Gland	M20 cable entry gland
Max. Temperature	150°C
Min. Temperature	-50°C
IP Rating	IP67
Transmitter	Optional transmitter device has an ambient operating range of -40 to +85°C
Number of Wires	4
Series	Outdoor / Cold Store Pt100