



Dtecteurs film plat 100 ohms - Classe A - 1,2 mm x 1,6 mm



Labfacility are currently multi-sourcing flat film detectors in order to ensure we can maintain a supply service to our customers without affecting the quality or functionality for the foreseeable future.

All flat film detectors supplied, although differing in appearance and packaging in some cases, will meet the specifications as published on our data sheets.



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

LabFacility fournit actuellement plusieurs fois des détecteurs film plat afin de garantir que nous pouvons maintenir un service d'approvisionnement nos clients sans affecter la qualité ou la fonctionnalité dans un avenir prévisible.

Tous les détecteurs film plat fournis, bien que différant dans certains cas par leur apparence et leur emballage, répondront aux spécifications publiées dans nos fiches techniques.

Caractéristiques techniques

Specifications

Product Code	XE-4440-001
Sensor Type	Pt100 (100 Ohms @ 0°C)
Element	Thin film, 10mm tails
Max. Temperature	500°C
Min. Temperature	-50°C
Tolerance	Class A and in accordance with IEC 60751
Dimensions	1.2mm Width x 1.6mm Length
Fundamental interval	38.5 ohms (nominal)
Ice point resistance	100 ohms
Self-heating	<0.5°C/mW
Stability	±0.05%
Thermal response	<0.4s @ t0.9 (in water at 0.4m/sec)
Important Information	Resists vibration & shock
Series	RTD Detectors