



Detecteurs film plat 100 Ohm - Classe 1/3 DIN - 2mm x 2.3mm



Labfacility are currently multi-sourcing flat film detectors in order to ensure we can maintain a supply service to our customers without affecting the quality or functionality for the foreseeable future.

All flat film detectors supplied, although differing in appearance and packaging in some cases, will meet the specifications as published on our data sheets.



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Labfacility sapprovisionne actuellement en dtecteurs couche plate afin de garantir que nous pouvons maintenir un service dapprovisionnement nos clients sans affecter la qualit ou la fonctionnalit dans un avenir prvisible.

Tous les dtecteurs film plat fournis, bien que diffrant dans certains cas en termes dapparence et demballage, rpondront aux spcifications publies sur nos fiches techniques.

Spcifications

Specifications

Product Code	XE-4466-001
Sensor Type	Pt100 (100 Ohms @ 0°C)
Element	Thin film, 10mm tails
Max. Temperature	500°C
Min. Temperature	-50°C
Tolerance	1/3DIN in accordance with IEC 60751
Dimensions	2mm Width x 2.3mm Length
Fundamental interval	38.5 ohms (nominal)
Ice point resistance	100 ohms
Self-heating	<0.5°C/mW
Stability	±0.05%
Thermal response	<0.4s @ t0.9 (in water at 0.4m/sec)
Important Information	Resists vibration & shock
Series	RTD Detectors