


Dtecteurs fil bobin 100 Ohm - Classe B - 1,5 mm x 15 mm



100 Ohm Fil Bobin 100 Ohm - Classe B - 1,5 mm x 15 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Veillez noter que les dlais de livraison DIN de classe 1/10 varient - veuillez vous renseigner ou vous serez contact avec un dlai aprs avoir pass votre commande.

Spcifications

Specifications

Product Code	XE-4710-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Pt100 Single Element
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	Class B
Dimensions	1.5mm Width x 15mm Length
Series	RTD Detectors