




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple bande magnétique avec fiche miniature Type K IEC - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple bande magnétique souple capable d'être fixée sur la plupart des surfaces ferromagnétiques, la grande surface de l'aimant assure une bonne adhérence permettant une simplicité d'installation. Utile dans un large éventail d'applications de surveillance de la température, en particulier dans les tests et le développement.

Specifications

Specifications

Product Code	XE-3455-001
General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Magnet	50 x 25 x 1mm magnetic strip
Cable Length	2M
Cable Type	PTFE insulated wires, twin twisted
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	100°C
Min. Temperature	-50°C
Standards Met	IEC
Series	Magnet Thermocouple