


**Thermocouple mis la terre fabriqué par CEI dans un tube en acier inoxydable, câble surtressé en acier inoxydable en fibre de verre - Type J - 4,76 mm de diamètre - 100 mm de longueur**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

|                              |                                              |
|------------------------------|----------------------------------------------|
| <b>Product Code</b>          | XE-5551-001                                  |
| <b>General Description</b>   | 316 stainless steel tube in various sizes    |
| <b>Thermocouple Type</b>     | J                                            |
| <b>Thermocouple Junction</b> | grounded at tip for fast response            |
| <b>Cable Length</b>          | 2M                                           |
| <b>Cable Type</b>            | Glassfibre stainless steel overbraided cable |
| <b>Cable Termination</b>     | Bare wire tails                              |
| <b>Max. Temperature</b>      | 350°C                                        |
| <b>Min. Temperature</b>      | -60°C                                        |
| <b>Standards Met</b>         | IEC                                          |