


## Thermomètre de résistance hygiénique Pt100, 1.5" RJT Style Fitting



[Click here to view product](#)

|                     |                                                                                        |
|---------------------|----------------------------------------------------------------------------------------|
| Product Code        | XE-5525                                                                                |
| General Description | Ideal for hygienic processes including dairy, brewery and pharmaceutical applications. |
| Sensor Type         | Pt100, 3 or 4 wire single element class A                                              |
| Head Type           | KNEF                                                                                   |
| Connection Type     | Pt100, 3 or 4 wire single element class A                                              |
| Transmitter         | Optional SEM-203p transmitter, ranged -50°C to +150°C                                  |
| Support Dia (D2)    | 8.0mm                                                                                  |
| Support Length (L2) | 50mm                                                                                   |