


**Thermocouple jonction expose CEI paire plate
isole en fibre de verre - Type K avec terminaison
de fiche miniature - 1/0,508 mm x 2 M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-5569-001
General Description	Good temperature resistance (not suitable in fluids)
Sensor Type	Exposed Junction
Thermocouple Type	K
Cable Length	2M
Cable Type	Glassfibre
Core / Strands	1/0.508mm
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
Tolerance	Class 2
Standards Met	IEC-584-3
Connector Type	Plug
Connector Size	Miniature
Colour	Green
Series	Exposed Junctions