



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple jonction expose CEI avec mini-prise intgre - Type K - PFA double torsad - Longueur 1 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple cble avec prise IEC miniature intgre, PFA double torsad ou paire plate.

Un thermocouple complet fil fin de prcision, jonction expose haute sensibilit, avec un connecteur miniature intgr.

Spcifications

Specifications

Product Code	XE-5598-001
General Description	A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector and cable tidy.
Sensor Type	Exposed Junction
Thermocouple Type	K IEC
Cable Length	1M
Cable Type	PFA Twin Twisted Insulated Wire
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	Rated 220°C continuous
Tolerance	Class 1
Standards Met	IEC 584
Connector Type	Plug
Connector Size	Miniature
Series	Cable Tidy Exposed Junction