

**Thermocouple jonction expose CEI avec  
mini-prise intgre - Type K - PFA double torsad -  
Longueur 2M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple cble avec prise IEC miniature intgre, PFA double torsad ou paire plate.

Un thermocouple complet fil fin de prcision, jonction expose haute sensibilit, avec un connecteur miniature intgr.

### **Vous avez besoin de calibrer votre capteur de temprature ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne, vous pouvez slectionner nimporte quel point de temprature entre -10C et +200C.

**Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison**

### **Spcifications**

## **Specifications**

<b>Product Code</b>	XE-7028-001
<b>General Description</b>	A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector and cable tidy.
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	K IEC
<b>Cable Length</b>	2M
<b>Cable Type</b>	PFA Twin Twisted Insulated Wire
<b>Termination Type</b>	Miniature flat pin plug, colour coded in accordance IEC 584
<b>Max. Temperature</b>	Rated 220°C continuous
<b>Tolerance</b>	Class 1
<b>Standards Met</b>	IEC 584
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Miniature