



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple jonction expose rang par cble ANSI PFA avec mini-prise intgre - Type K - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-5607-001
General Description	A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector
Sensor Type	Cable Tidy Thermocouple with Fitted Miniature ANSI Plug
Thermocouple Type	K
Thermocouple Junction	Exposed
Cable Length	2M
Cable Type	PFA twin twisted
Cable Termination	Miniature Plug
Max. Temperature	220°C
Tolerance	Class 1 wire for high accuracy
Standards Met	ANSI colour coded
Connector Type	Plug
Connector Size	Miniature
Colour	Yellow
Series	Cable Tidy Exposed Junction