

Thermocouple boulon IEC, acier inoxydable en fibre de verre sur cble tress - type J - pas de filetage M8x1mm - longueur de boulon de 13 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-5557-001
General Description	Bolt sensors are generally screwed into the surface which normally increases the speed of response and sensitivity of the measurement. Suitable for use up to +250°C
Thermocouple Type	J
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2M
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare tails
Max. Temperature	250°C
Standards Met	IEC
Bolt thread/pitch	M8 x 1mm