

## GM1311 Thermomtre alimentaire numrique professionnel haute performance



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**VEUILLEZ NOTER : Le dlai de livraison sera de 7 10 jours ouvrables lors du choix de l'etalonnage**

GM1311 Thermomtre numrique professionnel haute performance pour aliments/barbecue/viande

Temperature ambiante : 0 40 C (32 104 F)

Humidit de lenvironnement : 10-90% HR

Tension nominale : 3 V

1 x thermomtre alimentaire

1 x manuel d'utilisation en anglais

1 x CR2032

### Spfications

## Specifications

<b>Product Code</b>	XE-4002-001
<b>Accuracy</b>	-10 - 100°C: +/-1°C / 100 - 200°C: +/-1.5°C
<b>Resolution</b>	0.1°C / °F
<b>Colour</b>	White
<b>Dimensions</b>	9.29 in x 1.34 in x 0.71 in (23.6 cm x 3.4 cm x 1.8 cm). Display dimension:44x18mm
<b>Function #1</b>	LCD display
<b>Function #2</b>	Temperature lock
<b>Function #3</b>	Max/Min measuring
<b>Function #4</b>	Conversion between Celsius degree and Fahrenheit degree
<b>Function #5</b>	Auto power off
<b>Max. Temperature</b>	300°C (572°F)
<b>Min. Temperature</b>	-50°C (-58°F)
<b>Probe Length</b>	149.5 x 6mm



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

<b>Weight</b>	Weight: 0.95 oz (27 g)
<b>Power Supply</b>	CR2032 button battery x 1 (included)