


**Thermocouple aimant de traction IEC Type K de  
3 kg, cble isol PFA avec tresse en acier  
inoxydable se terminant en miniature ou en prise  
standard**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	Rapid, reliable temperature sensing on ferrous surfaces
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded measuring junction
<b>Magnet</b>	1" diameter (3kg pull)
<b>Cable Length</b>	1 or 2 metre
<b>Cable Type</b>	stainless steel over-braided PFA Teflon®
<b>Cable Termination</b>	Fitted Miniature or Standard Plug
<b>Max. Temperature</b>	+250°C (200°C Plug)