



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Thermocouple patch en silicone CEI, PFA avec prise miniature - Type K - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

|                            |                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | XE-5619-001                                                                                     |
| <b>General Description</b> | Resistant to many oils & chemicals                                                              |
| <b>Sensor Type</b>         | Moulded silicone patch thermocouple, 30 x 15mm patch, self adhesive with aluminium foil backing |
| <b>Thermocouple Type</b>   | K                                                                                               |
| <b>Cable Length</b>        | 2M                                                                                              |
| <b>Cable Type</b>          | PFA insulated twin-twist wire                                                                   |
| <b>Termination Type</b>    | Miniature flat pin plug, colour coded in accordance IEC 584                                     |
| <b>Max. Temperature</b>    | 150°C                                                                                           |
| <b>Min. Temperature</b>    | -50°C                                                                                           |
| <b>Standards Met</b>       | IEC                                                                                             |