

## **Thermocouple IEC expos double torsade isol en PTFE avec bouchon miniature - Type J - longueur 2M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple de type K, J ou T. Jonction apparente soude réponse rapide, disponible dans la spécification de câble de classe 1 conformément à la norme BS EN 60584-3:2008 avec une terminaison de prise miniature.

Ideal pour les applications de test et développement.

### Caractéristiques techniques

## Specifications

|                     |                                                            |
|---------------------|------------------------------------------------------------|
| Product Code        | XE-0419-001                                                |
| General Description | Fast response, welded exposed junction                     |
| Sensor Type         | Exposed Junction                                           |
| Thermocouple Type   | J IEC                                                      |
| Cable Length        | 2M                                                         |
| Cable Type          | Teflon® PTFE insulated twin-twisted lead                   |
| Core / Strands      | 1/0.2mm                                                    |
| Termination Type    | Miniature flat pin plug colour coded in accordance IEC 584 |
| Max. Temperature    | 260°C                                                      |
| Min. Temperature    | -75°C                                                      |
| Tolerance           | Class 1                                                    |
| Standards Met       | IEC-584-3                                                  |
| Connector Type      | Plug                                                       |
| Connector Size      | Miniature                                                  |