

Thermocouple double isolation minrale CEI avec prise duplex standard de type K - 6,0 mm x 500 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Vous avez besoin de calibrer votre capteur de température ?

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne.

Veuillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.

Spifications**Specifications**

Product Code	XE-2015-001
General Description	310 stainless steel sheath. Sheath can be bent/formed to suit many applications and processes.
Sensor Type	Mineral insulated Type 'K' with Duplex Thermocouple Plug
Thermocouple Type	K
Thermocouple Junction	Insulated hot junction
Probe Diameter	6.0mm
Probe Length	500mm
Termination Type	Standard duplex plug termination (200°C)
Max. Temperature	1100°C
Min. Temperature	-40°C
Standards Met	Conforms to IEC 584 specification
Connector Type	Plug
Connector Size	Standard