

**Sonde Pt100 usage gnral avec fil PFA - Type PRT**  
**- Diamtre de la sonde 3,0 mm - Longueur de la sonde 150 mm - Longueur du cble 2 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Sonde Pt100 usage gnral avec fil PFA - Type RTD / PRT**

Sonde polyvalente Pt100 selon IEC751

Classe B, connexion 4 fils

Gaine rigide en acier inoxydable 316

Plomb blind isol au tflon

Température nominale du capteur -50C +250C

**Spifications****Specifications**

|                            |                                                                       |
|----------------------------|-----------------------------------------------------------------------|
| <b>Product Code</b>        | XE-3609-001                                                           |
| <b>General Description</b> | Multipurpose Pt100 Probe to IEC751                                    |
| <b>Sensor Type</b>         | Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire                    |
| <b>Sheath Material</b>     | 3.0mm diameter 316 Stainless Steel Sheath                             |
| <b>Cable Length</b>        | 2M                                                                    |
| <b>Cable Type</b>          | Flexible Teflon® insulated, screened                                  |
| <b>Termination Type</b>    | 4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire |
| <b>Max. Temperature</b>    | 250°C                                                                 |
| <b>Min. Temperature</b>    | -50°C                                                                 |
| <b>Standards Met</b>       | IEC751                                                                |