


**Thermocouple IEC jonction expose paires isoles
isoles en PFA avec bouchon miniature - Type K -
Longueur du brin de 1/0,3 mm - 2 m de longueur**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bonne rsistance mcanique et flexibilit, rsistante aux huiles, acides et autres fluides indsirables.

Caractristiques techniques

Specifications

Product Code	XE-3485-001
Sensor Type	Fast response, welded exposed junction
Thermocouple Type	K IEC
Cable Connection	Terminated with a miniature flat pin plug
Cable Length	2M
Core / Strands	1/0.3mm
Termination Type	Miniature Plug
Standards Met	(cores & jacket) to IEC-584-3
Insulation Material	Insulation rating: -75°C to 260°C
Series	Exposed Junctions