


**Thermocouple jonction expose IEC paires isoles  
et isol en PFA avec bouchon miniature - Type K -  
Taille du brin de 1/0,3 mm - 3 M de longueur**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bonne rsistance mcanique et flexibilit, rsistante aux huiles, acides et autres fluides indsirables.

### Caractristiques techniques

## Specifications

<b>Product Code</b>	XE-3486-001
<b>Sensor Type</b>	Fast response, welded exposed junction
<b>Thermocouple Type</b>	K IEC
<b>Cable Connection</b>	Terminated with a miniature flat pin plug
<b>Cable Length</b>	3M
<b>Core / Strands</b>	1/0.3mm
<b>Termination Type</b>	Miniature Plug
<b>Standards Met</b>	(cores & jacket) to IEC-584-3
<b>Insulation Material</b>	Insulation rating: -75°C to 260°C
<b>Series</b>	Exposed Junctions