

Thermocouple jonction expose IEC PTFE avec prise miniature moule - Type K - Longueur 5M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un thermocouple isol en PTFE rponse trs rapide, fabriqu partir dune construction de cble unique qui fournit une isolation tanche au gaz, la vapeur et leau le long du cble.

Ce capteur est parfaitement adapt la surveillance de la temprature des pouttes dans les applications dautoclave et arospatiales.

Spifications

Specifications

Product Code	XE-0031-001
General Description	Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is wildly used in autoclaves.
Sensor Type	Exposed Junction
Thermocouple Type	K IEC
Cable Length	5M
Cable Type	PTFE
Core / Strands	1/0.376mm
Termination Type	Terminated with a miniature flat pin plug (Plug max temp 220°C)
Max. Temperature	250°C
Standards Met	IEC-584-3
Connector Type	Plug
Connector Size	Miniature