


**cran en cuivre plaqu tain 3 curs isol en PFA, cble  
/ fil de capteur PRT**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**cran en cuivre plaqu tain 3 curs isol en PFA, cble / fil de capteur PRT**

Temperature leve jusqu 260C, construction en tflon PFA extrud. Couramment utilis pour les thermomtres Pt100, galement adapt dautres types tels que Pt130, Pt500 et Pt1000

**Caractristiques techniques****Specifications**

<b>General Description</b>	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores, tin plated copper braided screen with PFA insulated outer jacket
<b>Cable Length</b>	10, 25, 50, 100m
<b>Cable Type</b>	PFA Teflon® construction
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	PFA Insulated
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-75°C
<b>Number of Wires</b>	3
<b>Standards Met</b>	Cores colour coded Red & White in accordance with IEC-751, black outer jacket