

Thermomètre rouge de sonde de pliage



[Click here to view product](#)

Product Code	XE-3951-001
Model	TP305
General Description	Household, folding thermometer, digital folding BBQ grill thermometer.
Type	Celsius and Fahrenheit
Accuracy	+/-1°C at -20°C to 200°C. +/-2°C at other ranges
Resolution	0.1 C/F
Response	Sampling time: 1 second / Response time: 4 seconds
Battery	AAA battery powered (included)
Dimensions	270mm x 35mm x 22mm. Probe length: 120mm / Probe diameter: 3mm
Functions	Fold / unfold the probe to switch on / off

Function #1	10 minutes auto-off
Function #2	Super fast response time
Probe Diameter	3mm
Probe Length	120mm
Probe Material	ABS and stainless steel
Sensor Type	Temperature Sensor
Sheath Material	ABS plastic. LCD digital display
Ranges Â°C	-50°C to 300°C