

**Thermocouple magnétique bouton de traction IEC  
de type K 0,7 kg, isol en PTFE avec prise  
standard de type K - Longueur de 1 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

|                              |                                             |
|------------------------------|---------------------------------------------|
| <b>Product Code</b>          | XE-3461-001                                 |
| <b>General Description</b>   | Ideal for ferrous metal surface measurement |
| <b>Thermocouple Type</b>     | K                                           |
| <b>Thermocouple Junction</b> | Grounded at tip                             |
| <b>Magnet</b>                | 1/2" diameter button type                   |
| <b>Cable Length</b>          | 1M                                          |
| <b>Cable Type</b>            | PTFE insulated twin-twist wires             |
| <b>Cable Termination</b>     | Standard Fitted Plug                        |
| <b>Max. Temperature</b>      | +250°C (plug 200°C)                         |