



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple isolation minrale de style industriel avec tte de raccord KNE - Inconel de type K (6 mm x 400 mm)

Thermocouple isolation minrale de style industriel avec tte de raccord KNE - Inconel de type K (6 mm x 400 mm)



6.0 X 100 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-3248-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral insulated Thermocouple fitted with KNE terminal head
Thermocouple Type	Type K
Thermocouple Junction	Insulated hot junction
Probe Diameter	6.0mm diameter
Probe Length	400mm Length
Sheath Material	Inconel
Head Type	KNE, Cast aluminium finish, screw type lid with chain
Gland	M16 x 1.5mm Plated brass gland fitted
Max. Temperature	+1100°C
Min. Temperature	-40°C
IP Rating	IP67