



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Thermocouple isolation minrale rponse rapide
ANSI de 0,25 mm de diamtre avec prise miniature
ANSI - Jonction mise la terre de type K -
Longueur de 100 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

[formulaires de capteurs personnalisés](#)

Veuillez utiliser nos
pour soumettre vos besoins pour différentes longueurs de MI et de longueur

Vous avez besoin de calibrer votre capteur de température ?

Traable selon les normes d'étalonnage ISO17025 (UKAS), Labfacility offre un choix d'étalonnage traable en 3 ou 5 points laide de notre installation d'étalonnage interne.

Veuillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.

Spécifications

Specifications

Product Code	XE-0809-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Type K ANSI MB STYLE probes 0.25mm
Element/hot junction	Single element, grounded junction
Probe Length	100mm
Sheath Material	Stainless steel 316 sheath
Termination Type	Fitted Miniature Plug Type ANSI
Max. Temperature	260°C
Min. Temperature	-40°C
Tolerance	A 1%+ Length Tolerance applies
Standards Met	ANSI MC96.1
Connector Type	Plug
Connector Size	Miniature



Colour	Yellow
--------	--------