

0.25mm Diameter Fast Response Mineral Insulated Thermocouple with ANSI Miniature Plug - Type K



IMAGE
COMING
SOON

[Click here to view product](#)

Product Code	XE-0809
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Type K ANSI MB STYLE probes 0.25mm
Thermocouple Junction	Single element, junction insulated or grounded
Probe Length	100mm / 150mm / 250mm
Sheath Material	Stainless steel 316 sheath
Termination Type	Fitted Miniature Plug Type ANSI
Range	-40 to +260°C
Important Information	A 1%+ Length Tolerance applies