



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Thermocouple isolation minrale rponse rapide
de 0,25 mm de diamtre avec pot filet et queues
IEC PFA - Type K - Jonction isole - 100 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Vous avez besoin de calibrer votre capteur de température ?

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne, vous pouvez slectionner nimporte quel point de temprature entre -10C et +200C.

Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.

Specifications

| | |
|------------------------------|--|
| Product Code | XE-0912-001 |
| General Description | Highly flexible, sheath can be bent/formed to suit many applications and processes |
| Sensor Type | Type K IEC MD STYLE probes 0.25mm |
| Thermocouple Junction | Single element, junction insulated |
| Probe Length | 100mm |
| Sheath Material | Stainless steel 316 sheath |
| Pot Type | M8 x 1.0mm fine pitch threaded pot seal (200°C) |
| Cable Length | 100mm |
| Cable Type | PFA insulated twin twisted IEC |
| Core / Strands | 7/0.2mm |
| Range | Sensor range -40 to +260°C |
| Important Information | A 1%+ Length Tolerance applies |