

**Thermocouple isolation minrale ANSI de 0,5 mm
de diamtre avec pot et queues filets ANSI -
Jonction isole - Longueur de la sonde de 100 mm**

ANSI 100mm MINERAL INSULATED THERMOCOUPLE
WITH POT AND ANSI WIRE LEADS



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Type K avec pot fileté et queues en acier inoxydable 310

Pour les diamètres de 0,25 mm, 0,75 mm et 0,81 mm et pour les versions de gaine en alliage Inconel 600, veuillez utiliser nos

[formulaires de capteurs personnalisés](#)

pour soumettre vos besoins ou envoyer un e-mail sales@labfacility.com / téléphone notre 871280

Vous avez besoin de calibrer votre capteur de température ?

Traçable selon les normes d'étalonnage ISO17025 (UKAS), Labfacility offre un choix d'étalonnage traçable en 3 ou 5 points. À l'aide de notre installation d'étalonnage interne, vous pouvez sélectionner n'importe quel point de température entre -10°C et +200°C.

Veuillez noter que si vous ajoutez une option de calibrage à votre commande, vous serez contacté pour les informations de livraison.

Spécifications

Specifications

Product Code	XE-0181-001
General Description	Type K with threaded pot & tails in 310 stainless steel
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction insulated from sheath
Probe Diameter	0.5mm
Probe Length	100mm
Sheath Material	310 stainless steel
Pot Type	M8 x 1.0mm fine pitch threaded pot seal (200°C)
Cable Length	100mm
Cable Type	PFA insulated twin twisted ANSI
Core / Strands	7/0.2mm



Termination Type	Bare tails
Max. Temperature	+650°C
Min. Temperature	-40°C
Standards Met	In accordance with ANSI MC96.1