



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Thermocouple réponse rapide de 0,75 mm de
diamètre avec pot filet et queues IEC - Jonction
isole - 250 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Type K avec pot et queues filets en acier inoxydable 316

Si vous avez besoin de plus grandes quantités, veuillez sales@labfacility.com envoyer un e-mail ou téléphoner notre ligne de vente 01243 871280.

Vous avez besoin de calibrer votre capteur de température ?

Traçable selon les normes d'étalonnage ISO17025 (UKAS), Labfacility offre un choix d'étalonnage traçable en 3 ou 5 points. À l'aide de notre installation d'étalonnage interne, vous pouvez sélectionner n'importe quel point de température entre -10°C et +200°C.

Veuillez noter que si vous ajoutez une option de calibrage à votre commande, vous serez contacté pour les informations de livraison.

Spécifications

Specifications

Product Code	XE-0202-001
General Description	Type K with threaded pot & tails in either 316 stainless steel.
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction insulated from sheath
Probe Diameter	0.75mm
Probe Length	250mm
Sheath Material	316 stainless steel
Pot Type	M8 x 1.0mm fine pitch threaded pot seal (Max Temp 200°C)
Cable Length	100mm
Cable Type	PFA insulated twin twisted IEC (Temp -75 to +250°C)
Core / Strands	7/0.2mm
Termination Type	Bare Tails



Max. Temperature	+750°C
Min. Temperature	-100°C
Tolerance	Class 1
Standards Met	In accordance with IEC 584