



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple réponse rapide de 0,75 mm de diamètre avec queues nues IEC - Jonction isolée - 100 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sondes thermocouples joint de transition de diamtre fin

Les jonctions de thermocouple sont disponibles en version mise la terre, pour un temps de rponse plus rapide et un cot infrieur ; ou isol, pour lisolation lectrique du thermocouple de la gaine.

Si vous avez besoin de plus grandes quantits, veuillez sales@labfacility.com envoyer un e-mail ou tlphoner notre ligne de vente 01243 871280.

Vous avez besoin de calibrer votre capteur de temprature ?

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne, vous pouvez slectionner nimporte quel point de temprature entre -10C et +200C.

Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.

Spfications**Specifications**

Product Code	XE-0203-001
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction insulated from sheath
Probe Diameter	0.75mm
Probe Length	100mm
Sheath Material	K - Stainless steel 316
Termination Type	Fortfix sealed with = +VE leg 15mm long and -VE leg 10mm long
Max. Temperature	+650°C
Min. Temperature	-40°C
Tolerance	Class 1



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Standards Met	In accordance with IEC 584
----------------------	----------------------------