



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple réponse rapide de 0,75 mm de diamtre avec queues nues IEC - Jonction isole - 250 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sondes thermocouples joint de transition de diamètre fin

Les jonctions de thermocouple sont disponibles en version mise à la terre, pour un temps de réponse plus rapide et un coût inférieur ; ou isolé, pour l'isolation électrique du thermocouple de la gaine.

Si vous avez besoin de plus grandes quantités, veuillez sales@labfacility.com envoyer un e-mail ou téléphoner notre ligne de vente 01243 871280.

Vous avez besoin de calibrer votre capteur de température ?

Traçable selon les normes d'étalonnage ISO17025 (UKAS), Labfacility offre un choix d'étalonnage traçable en 3 ou 5 points. À l'aide de notre installation d'étalonnage interne, vous pouvez sélectionner n'importe quel point de température entre -10°C et +200°C.

Veuillez noter que si vous ajoutez une option de calibrage à votre commande, vous serez contacté pour les informations de livraison.

Spécifications

Specifications

Product Code	XE-0205-001
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction insulated from sheath
Probe Diameter	0.75mm
Probe Length	250mm
Sheath Material	K - Stainless steel 316
Termination Type	Fortifix sealed with = +VE leg 15mm long and -VE leg 10mm long
Max. Temperature	+650°C
Min. Temperature	-40°C
Tolerance	Class 1



Standards Met

In accordance with IEC 584