



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple minral isol minral de 0,81 mm de diamtre avec fichire ANSI miniature jonction mise la terre - 250 mm de longueur de la sonde



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Type K avec prise ANSI miniature en acier inoxydable 316

Si vous avez besoin de plus grandes quantités, veuillez sales@labfacility.com envoyer un courriel ou appeler notre ligne commerciale au 01243 871280.

Vous avez besoin d'étalonner votre capteur de température ?

Traable selon les normes d'étalonnage ISO17025 (UKAS), Labfacility propose le choix entre un étalonnage traable en 3 ou 5 points via notre installation d'étalonnage interne.

Veuillez noter que si vous ajoutez une option de calibration votre commande, vous serez contacté pour les informations de livraison.

Specifications

Specifications

Product Code	XE-0253-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction grounded from sheath
Probe Diameter	0.81mm
Probe Length	250mm
Sheath Material	K - Stainless steel 316
Pot Type	Plain (+200°C)
Cable Length	1M
Cable Type	PFA Teflon® insulated flat pair cable and tails (colour coded to ANSI MC96.1) (Max Temp 250°C)
Core / Strands	7/0.2mm



Termination Type	Miniature flat pin plug, colour coded In accordance with ANSI MC96.1
Max. Temperature	+650°C. Plug Max. temperature range: 200°C
Min. Temperature	-40°C
Tolerance	Class 1
Standards Met	ANSI MC96.1
Connector Type	Plug
Connection Type	Miniature
Colour	Yellow