


**Thermocouple magnétique IEC Type K 11,8 kg (fer cheval), cble isol PFA avec sur-tresse en acier inoxydable se terminant par une prise standard - 3,0 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Spcifications****Specifications**

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Product Code</b>        | XE-0078-001                                         |
| <b>General Description</b> | 11.8kg Pull Power (Horseshoe) Magnet Thermocouple   |
| <b>Sensor Type</b>         | Magnet Thermocouple                                 |
| <b>Thermocouple Type</b>   | K                                                   |
| <b>Magnet</b>              | 45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG      |
| <b>Cable Length</b>        | 3.0M                                                |
| <b>Cable Type</b>          | PFA Insulated Cable with Stainless Steel Over-Braid |
| <b>Core / Strands</b>      | 7/0.2mm                                             |
| <b>Termination Type</b>    | Standard Plug                                       |
| <b>Max. Temperature</b>    | 250°C                                               |
| <b>Min. Temperature</b>    | -50°C                                               |
| <b>Standards Met</b>       | IEC                                                 |