


Connecteur thermocouple standard Douille en ligne FSTC-CU-F Type Cu



Standard in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur thermocouple standard Douille en ligne FSTC-CU-F Type CU

Accepte les broches rondes standard, adaptes aux applications industrielles

Caractristiques techniques**Specifications**

Product Code	XE-1672-001
General Description	Accepts standard round plug pins, contacts are polarised to prevent incorrect connection.
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Socket
Connector Size	Standard
Colour	White
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS