

Thermocouple boulon ANSI, acier inoxydable en fibre de verre sur cble tress - type K - filetage M8 - boulon de 25 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0585-001
Thermocouple Type	K
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2M
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare tails
Max. Temperature	250°C
Standards Met	ANSI MC96.1
Bolt thread/pitch	M8 x 1mm
Colour	Positive: Red / Negative: Yellow
Length	25mm Bolt
Series	Bolt Thermocouple